

Ротаторы MX-RD-E, MX-RL-E, шейкеры для обесцвечивания SK-L180-S, SK-L180-Pro, SK-R1807-S, 3D-шейкеры SK-D1810-S

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Ноябрьск (3496)41-32-12
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: dbe@nt-rt.ru || сайт: <http://dlab.nt-rt.ru/>

MX-RD-E Effective affordable mixing

- Speed controlled by Scale
- Adjustable speed range from 0 to 80 rpm

Model : **MX-RD-E**

Cat.No.

8031412000

Accessories type

Disk

Tilt angle

-

Motion type

-

Rotation angle

-

Amplitude

-

Speed range

0-80rpm

Speed display

Scale

Timer display

-

Time setting range

-

Run type

Continuous

Max. load capacity

-

Number of rollers

-

Roller size[length]

-

Motor type

DC motor

Voltage

100-240V

Frequency

50/60Hz

Power

20W

Dimension [D×W×H]

300×220×310mm

Weight

-

Permissible ambient temperature

5-40°C

Permissible ambient humidity

80%RH

Protection class

IP21

MX-RL-E Economical Rotator

- Provides gentle but effective rolling
- Adjustable speed range from 0 to 80 rpm
- Used with MX-RL-E/MX-RD-E Disk/Rotisserie accessory, for 50mL, 15mL and 1.5mL

Model : **MX-RL-E**

PRODUCT PARAMETERS

Cat.No.	8031402000
Accessories type	Tube Rotator
Tilt angle	-
Motion type	-
Rotation angle	-
Amplitude	-
Speed range	0-80rpm
Speed display	Scale
Timer display	-
Time setting range	-
Run type	Continuous
Max. load capacity	-
Number of rollers	-
Roller size[length]	-
Motor type	DC motor
Voltage	100-240V
Frequency	50/60Hz
Power	20W
Dimension [D×W×H]	150X530X190mm
Weight	-
Permissible ambient temperature	5-40°C
Permissible ambient humidity	80% RH
Protection class	

SK-L180-S Intelligent Linear Decolorizing Shaker

- Ensure uniform mixing to enhance experiment reliability and accuracy.
- Provide different load capacities suitable for experiments of varying scales and sample processing, meeting diverse experimental needs.
- Independent display screen, simultaneously showing speed and time, user-friendly operation for easy control of the experiment process.
- Offer a wide range of speed settings that can be adjusted according to specific experimental requirements to ensure optimal results.
- Maintenance-free brushless DC motor for extended lifespan.
- Equipped with safety speed protection features.

Suitable for mixing in culture dishes, flasks and beakers. Ideal for aeration and extraction applications. It can be used in environmental chambers and CO2 incubators.

Model : SK- L180-S

PRODUCT PARAMETERS

Cat.No.	8032130200
Motion type	Linear Motion
Amplitude	20mm
Tilt angle	-
Speed range	40-200rpm
Speed display	LED
Timer display	LED
Time settings range	1min-99h59min
Data connector	-
Max. load capacity (with platform)	2kg
Motor type	Brushless DC motor
Motor rating input	16W
Motor rating output	10W
Power	20W
Voltage	100-240V
Frequency	50/60Hz
Dimension[W×D×H]	294X348X123mm
Work plate Dimension[W x D]	268x268mm
Weight	3.4kg
Permissible ambient temperature	5-40°C
Permissible ambient humidity	80% RH
Protection class	

SK-L180-Pro Multifunctional Linear Decolorizing Shaker

- Linear reciprocating shaking action
- 2.5kg with platform
- Independent LCD display for speed and time
- Speed range of 100-350rpm
- Maintenance-free brushless DC motor
- Over speed detection and protection
- Wide range of platforms for use with a variety of vessels

Model : **SK-L180-Pro**

PRODUCT PARAMETERS

Cat.No.	8032110100
Motion type	Linear Shaker
Amplitude	4mm
Tilt angle	-
Speed range	100-350rpm
Speed display	LCD
Timer display	LCD
Time settings range	1min-19h59min
Data connector	RS232
Max. load capacity (with platform)	2.5kg
Motor type	Brushless DC motor
Motor rating input	28W
Motor rating output	15W
Power	30W
Voltage	100-240V
Frequency	50/60Hz
Dimension[W×D×H]	300X340X100mm
Work plate Dimension[W x D]	250×245mm
Weight	8.1kg
Permissible ambient temperature	5-40°C
Permissible ambient humidity	80% RH
Protection class	IP21

SK-R1807-S Intelligent Transference Decolorizing Shaker

- Rocking motion angle 7°=
- 2kg and speed range 10-80rpm
- Double LEDs for independent displays of time and speed
- Continuous or timer operation mode
- Electronic time switching clock: 1min-99h59min
- Fixed work plate
- Simple appearance, save operating space

Model : SK-R1807-S

PRODUCT PARAMETERS

Cat.No.	8072311200
Motion type	Rocking Shaker
Amplitude	-
Tilt angle	7°
Speed range	10-80rpm
Speed display	LED
Timer display	LED
Time settings range	1min-99h59min
Data connector	-
Max. load capacity (with platform)	2kg
Motor type	DC motor
Motor rating input	20W
Motor rating output	16W
Power	20W
Voltage	100-240V
Frequency	50/60Hz
Dimension[W×D×H]	294×348×123mm
Work plate Dimension[W x D]	268×268mm
Weight	3.1kg
Permissible ambient temperature	5-40°C
Permissible ambient humidity	80% RH
Protection class	

SK-D1810-S Intelligent 3D Shaker

- Provide a gentle three-dimensional motion shaker at a fixed tilt angle
- Anti slip silicone disc surface, better anti slip effect, and better sample processing effect
- Adjustable rotational mixing speed
- Independent LCD/LED display of time and speed
- Incubator compatible

3D shaker used for uniform shaking mixing in flasks, culture flasks, petri dishes and tubes. Ideal for culturing a variety of cell cultures, DNA extractions, distributing of cells and gel staining & destaining of cells.

Model : **SK-D1810-S**

PRODUCT PARAMETERS

Cat.No.	8032434200
Motion type	3D Motion
Amplitude	-
Tilt angle	10°
Speed range	10-80rpm
Speed display	LED
Timer display	LED
Time settings range	1min-99h59min
Data connector	-
Max. load capacity (with platform)	3kg
Motor type	Brushless DC motor
Motor rating input	16W
Motor rating output	10W
Power	20W
Voltage	100-240V
Frequency	50/60Hz
Dimension[W×D×H]	294X348X123mm
Work plate Dimension[W x D]	268×268mm
Weight	3.7kg
Permissible ambient temperature	5-40°C
Permissible ambient humidity	80% RH
Protection class	IP21

SK-L180-Pro

Multifunctional Linear Decolorizing Shaker

SK-L180-S

Intelligent Linear Decolorizing Shaker

Features

- Ensure uniform mixing to enhance experiment reliability and accuracy.
- Provide different load capacities suitable for experiments of varying scales and sample processing, meeting diverse experimental needs.
- Independent display screen, simultaneously showing speed and time, user-friendly operation for easy control of the experiment process.
- Offer a wide range of speed settings that can be adjusted according to specific experimental requirements to ensure optimal results.
- Maintenance-free brushless DC motor for extended lifespan.
- Equipped with safety speed protection features.
- Provide various platform options to accommodate a variety of container types and sizes, increasing flexibility and versatility (SK-L180-Pro).



SK-L180-Pro



SK-L180-S




Linear Motion

Applications

Suitable for mixing in culture dishes, flasks and beakers. Ideal for aeration and extraction applications. It can be used in environmental chambers and CO2 incubators.

Specifications	SK-L180-Pro	SK-L180-S
Amplitude	4mm	20mm
Max. load capacity(with platform)	2.5kg	2kg
Motor type	Brushless DC motor	Brushless DC motor
Speed range	100-350rpm	40-200rpm
Speed display	LCD	LED
Timer display	LCD	LED
Time settings range	1min-99h59min	1min-99h59min
Voltage, Frequency	100-240V, 50/60Hz	100~240V, 50/60Hz
Power	30W	20W
Motor input power	28W	16W
Motor output power	15W	10W
Weight	8.1kg	3.4kg
Permissible ambient temperature and humidity	5-40°C, 80%RH	5-40°C, 80%RH
Dimension[W×D×H]	300X340X100mm	294X348X130mm
Protection class	IP21	IP20
Data connector	RS232	-

Cat. No.	Description	SK-0300-Pro	SK-L330-Pro	SK-R330-Pro	SK-D3309-Pro	SK-0330-M	SK-01880-Pro	SK-L180-Pro	SK-0180-S (Multifunc)	SK-R180-S (Multifunc)
 18900155	SK330.5, Dish platform with non-slip mat and fixed string (8pcs) Dimension[WxD]: 30x28.5cm	■	■	■	■					
 18202835	Silicone mattress for 16× 5~15mL tubes, use with SK330.5	■	■	■	■					
 18202846	Silicone mattress for 6× 50mL tubes, use with SK330.5	■	■	■	■					
 18900185	SK180.4, Dish attachment with non-slip mat Dimension[WxD]: 24x21.5cm						■	■	■	■
 18900869	SK330.4 Dish platform with non-slip mat and fixed string(8 pcs), 34X34cm					■				
 18900859	SK330.6 Transparent shield, serologic test for syphilis, 41X39.5X10cm	■	■			■				
 18900860	SK180.5 Transparent Cover for Syphilis Serum Test, 32X32X8.3cm						■	■		
 18202852	Silicone mattress for 10× 5~10mL tubes						■	■	■	■
 18202853	Silicone mattress for 4× 50mL tubes						■	■	■	■
 18900027	SK330.1, Universal attachment with 4 bars Size: 32x32cm Capacity for flask: 25mLx21pcs, 50mLx15pcs, 100mLx12pcs, 250mLx6pcs, 500x4pcs	■	■	■	■					
 18900040	SK330.3 Lengthways roller attachment, Mainly used for separation funnel oscillation extraction applications 500mL*3pcs 250mL*4pcs 100/125mL*4pcs 50mL*5pcs	■	■	■	■					
 18900025	SK180.1, Universal attachment with 3 bars Size: 25x24.5cm Capacity for flask: 25mLx10pcs, 50mLx8pcs, 100mLx6pcs, 250mLx4pcs, 500x2pcs						■	■	■	■
 18900038	SK180.3 Lengthways roller attachment, Mainly used for separation funnel oscillation extraction applications 500mL*2pcs 250mL*3pcs 100/125mL*3pcs 50mL*4pcs						■	■	■	■

Cat. No.	Description	SK-O330-Pro	SK-L330-Pro	SK-R330-Pro	SK-D3309-Pro	SK-O330-M	SK-O180-Pro	SK-L180-Pro	SK-O180-S (Multifunc)	SK-R1807-S (Multifunc)
 18900036	Bar for SK330.1 universal attachment	■	■	■	■					
 18900037	Bar for SK180.1 universal attachment						■	■	■	■
 18900280	Connect sticks for double dishes shaker	■	■	■	■	■	■	■	■	■
 18900028	SK330.2, Fixing clip attachment Size: 32 x 32cm	■	■	■	■					
 18900026	SK180.2, Fixing clip attachment Size: 25 x 24cm						■	■	■	■
 18900029	Fixing clip for flask volume 25mL, used with SK330.2[28pcs] SK180.2[9pcs]	■	■	■	■		■	■	■	■
 18900030	Fixing clip for flask volume 50mL, used with SK330.2[16pcs] SK180.2[9pcs]	■	■	■	■		■	■	■	■
 18900031	Fixing clip for flask volume 100mL, used with SK330.2[16pcs] SK180.2 [9pcs]	■	■	■	■		■	■	■	■
 18900032	Fixing clip for flask volume 200/250mL, used with SK330.2[9pcs] SK180.2[4pcs]	■	■	■	■		■	■	■	■
 18900033	Fixing clip for flask volume 500mL, used with SK330.2[9pcs] SK180.2[2pcs]	■	■	■	■		■	■	■	■

SK-8330-Pro

Large Multifunctional Transference Rocking Shaker

SK-81801-S

Intelligent Transference Decolorizing Shaker

Features

- Rocking motion angle 9° (R330-Pro) and 7° (R1807-S)
- Two options:
10kg (R330-Pro) and speed range 10-70rpm,
2kg (R1807-S) and speed range 10-50rpm
- Double LCDs/LEDs for independent display of time
and speed
- Continuous or timer operation mode
- Electronic time switching clock: 1min-99h59min
- Wide range of platforms for use with a variety of
vessels(R330-Pro)
- Fixed work plate(R1807-S)
- Simple appearance, save operating space

Applications

Designed for time-dependent and low rpm sensitive shaking applications such as ELISA washes, enzyme immunoassays, protein synthesis, hybridization staining/destaining gels, immune precipitation and Western blots, cell culture, pharmaceutical profiling, etc.



Rocking Motion




SK-8330-Pro



SK-81807-S



SK-81807-S fixed plate



SK-81807-S multiple plate are optional

Specifications	SK-R330-Pro	SK-R1807-S
Tilt angle	9°	7°
Max. load capacity(with platform)	10kg	2kg
Motor type	DC motor	DC motor
Motor input power	40W	20W
Motor output power	24W	15W
Speed range	10-70rpm	10-80rpm
Speed display	LCD	LCD
Timer display	LCD	LCD
Time settings range	1min-99h59min	1min-99h59min
Voltage,Frequency	100-240V, 50/60Hz	100-240V, 50/60Hz
Power	50W	20W
Weight	9kg	3.1kg
Permissible ambient temperature and humidity	5-40°C, 80% RH	5-40°C, 80% RH
Dimension[W X D x H]	360 x 410 x 200mm	294 X 348 X 123mm
Protection class	IP21	IP21

Cat. No.	Description	SK-0330-Pro	SK-L330-Pro	SK-R330-Pro	SK-D3309-Pro	SK-0330-M	SK-0180-Pro	SK-L180-Pro	SK-0180-S (Multifunc)	SK-R807-S (Multifunc)
18900155	SK330.5, Dish platform with non-slip mat and fixed string (8pcs) Dimension[WxD]: 30x28.5cm	■	■	■	■					
18202835	Silicone mattress for 16X 5~15mL tubes, use with SK330.5	■	■	■	■					
18202846	Silicone mattress for 6X 50mL tubes, use with SK330.5	■	■	■	■					
18900185	SK180.4, Dish attachment with non-slip mat Dimension[WxD]: 24x21.5cm						■	■	■	■
18900869	Dish platform with non-slip mat and fixed string(8 pcs), 34X34cm					■				
18900859	SK330.6 Transparent shield, serologic test for syphilis, 41X39.5X10cm	■	■			■				
18900860	SK180.5 Transparent Cover for Syphilis Serum Test, 32X32X8.3cm						■	■		
18202852	Silicone mattress for 10X 5~10mL tubes						■	■	■	■
18202853	Silicone mattress for 4x 50mL tubes						■	■	■	■
18900027	SK330.1, Universal attachment with 4 bars Size: 32x32cm Capacity for flask: 25mLx21pcs, 50mLx15pcs, 100mLx12pcs, 250mLx6pcs, 500x4pcs	■	■	■	■					
18900040	SK330.3 Lengthways roller attachment, Mainly used for separation funnel oscillation extraction applications 500mL*3pcs 250mL*4pcs 100/125mL*4pcs 50mL*5pcs	■	■	■	■					
18900025	SK180.1, Universal attachment with 3 bars Size: 25x24.5cm Capacity for flask: 25mLx10pcs, 50mLx8pcs, 100mLx6pcs, 250mLx4pcs, 500x2pcs						■	■	■	■
18900038	SK180.3 Lengthways roller attachment, Mainly used for separation funnel oscillation extraction applications 500mL*2pcs 250mL*3pcs 100/125mL*3pcs 50mL*4pcs						■	■	■	■

Cat. No.	Description	SK-O330-Pro	SK-I330-Pro	SK-R330-Pro	SK-D3309-Pro	SK-O330-M	SK-O180-Pro	SK-L180-Pro	SK-O180-S (Multilunc)	SK-R1807-SI (Multilunc)
18900036	Bar for SK330.1 universal attachment	■	■	■	■					
18900037	Bar for SK180.1 universal attachment						■	■	■	■
18900280	Connect sticks for double dishes shaker	■	■	■	■	■	■	■	■	■
18900028	SK330.2, Fixing clip attachment Size: 32 x 32cm	■	■	■	■					
18900026	SK180.2, Fixing clip attachment Size: 25 x 24cm						■	■	■	■
18900029	Fixing clip for flask volume 25mL, used with SK330.2[28pcs] SK180.2[9pcs]	■	■	■	■	■	■	■	■	■
18900030	Fixing clip for flask volume 50mL, used with SK330.2[16pcs] SK180.2[9pcs]	■	■	■	■	■	■	■	■	■
18900031	Fixing clip for flask volume 100mL, used with SK330.2[16pcs] SK180.2[9pcs]	■	■	■	■	■	■	■	■	■
18900032	Fixing clip for flask volume 200/250mL, used with SK330.2[9pcs] SK180.2[4pcs]	■	■	■	■		■	■	■	■
18900033	Fixing clip for flask volume 500mL, used with SK330.2[9pcs] SK180.2[2pcs]	■	■	■	■		■	■	■	■

SK-D3309-Pro

Large Multifunctional 3D Shaker

SK-D1810-S

Intelligent 3D Shaker

Features

- Provide a gentle three-dimensional motion shaker at a fixed tilt angle
- Anti slip silicone disc surface, better anti slip effect, and better sample processing effect
- Adjustable rotational mixing speed
- Independent LCD/LED display of time and speed
- Incubator compatible
- Various platform accessories (D3309-Pro)

Applications

3D shaker used for uniform shaking mixing in flasks, culture flasks, petri dishes and tubes. Ideal for culturing a variety of cell cultures, DNA extractions, distributing of cells and gel staining & destaining of cells.



Wide range of accessories



3D Motion



SK-D3309-Pro



SK-D1810-S

Specifications	SK-D3309-Pro	SK-D1810-S
Tilt angle	9°	10°
Max. load capacity(with platform)	5kg	2kg
Motor type	DC motor	DC motor
Motor input power	40W	16W
Motor output power	24W	10W
Speed range	10-70rpm	10-50rpm
Speed display	LCD	LED
Timer display	LCD	LED
Time settings range	1min-99h59min	1min-99h59min
Voltage,Frequency	100-240V, 50/60Hz	100-240V,50/60Hz
Power	50W	20W
Weight	9kg	3.7kg
Permissible ambient temperature and humidity	5-40 °C,80%RH	5-40 °C,80%RH
Dimension[W x D x H]	360X430X 106mm	294X348X1 70mm
Protection class	IP21	IP21



Cat. No.	Description	SK-O330-Pro	SK-L330-Pro	SK-R330-Pro	SK-D3309-Pro	SK-O330-M	SK-O180-Pro	SK-L180-Pro	SK-O180-S (Multifunc)	SK-R1807-S (Multifunc)
 18900155	SK330.5, Dish platform with non-slip mat and fixed string (8pcs) Dimension[WxD]: 30x28.5cm	■	■	■	■					
 18202835	Silicone mattress for 16X 5~15mL tubes, use with SK330.5	■	■	■	■					
 18202846	Silicone mattress for 6X 50mL tubes, use with SK330.5	■	■	■	■					
 18900185	SK180.4, Dish attachment with non-slip mat Dimension[WxD]: 24x21.5cm						■	■	■	■
 18900869	Dish platform with non-slip mat and fixed string(8 pcs), 34X34cm					■				
 18900859	SK330.6 Transparent shield, serologic test for syphilis, 41X39.5X10cm	■	■			■				
 18900860	SK180.5 Transparent Cover for Syphilis Serum Test, 32X32X8.3cm						■	■		
 18202852	Silicone mattress for 10X 5~10mL tubes						■	■	■	■
 18202853	Silicone mattress for 4x 50mL tubes						■	■	■	■
 18900027	SK330.1, Universal attachment with 4 bars Size: 32x32cm Capacity for flask: 25mLx2(pcs, 50mLx3(pcs, 100mLx12(pcs, 250mLx6(pcs, 500x4(pcs	■	■	■	■					
 18900040	SK330.3 Lengthways roller attachment, Mainly used for separation funnel oscillation extraction applications 500mL*3pcs 250mL*4pcs 100/125mL*4pcs 50mL*5pcs	■	■	■	■					
 18900025	SK180.1, Universal attachment with 3 bars Size: 25x24.5cm Capacity for flask: 25mLx10(pcs, 50mLx8(pcs, 100mLx6(pcs, 250mLx4(pcs, 500x2(pcs						■	■	■	■
 18900038	SK180.3 Lengthways roller attachment, Mainly used for separation funnel oscillation extraction applications 500mL*2pcs 250mL*3pcs 100/125mL*3pcs 50mL*4pcs						■	■	■	■

Cat. No.	Description	SK-O330-Pro	SK-L330-Pro	SK-R330-Pro	SK-D3309-Pro	SK-O330-M	SK-O180-Pro	SK-L180-Pro	SK-O180-S (Multifunc)	SK-R1807-S (Multifunc)
 18900036	Bar for SK330.1 universal attachment	■	■	■	■					
 18900037	Bar for SK180.1 universal attachment						■	■	■	■
 18900280	Connect sticks for double dishes shaker	■	■	■	■	■	■	■	■	■
 18900028	SK330.2, Fixing clip attachment Size: 32 x 32cm	■	■	■	■					
 18900026	SK180.2, Fixing clip attachment Size: 25 x 24cm						■	■	■	■
 18900029	Fixing clip for flask volume 25mL, used with SK330.2(28pcs) SK180.2(9pcs]	■	■	■	■		■	■	■	■
 18900030	Fixing clip for flask volume 50mL, used with SK330.2(16pcs) SK180.2(9pcs]	■	■	■	■		■	■	■	■
 18900031	Fixing clip for flask volume 100mL, used with SK330.2(16pcs) SK180.2(9pcs]	■	■	■	■		■	■	■	■
 18900032	Fixing clip for flask volume 200/250mL, used with SK330.2(9pcs) SK180.2(4pcs]	■	■	■	■		■	■	■	■
 18900033	Fixing clip for flask volume 500mL, used with SK330.2(9pcs) SK180.2(2pcs]	■	■	■	■		■	■	■	■

MX-RD-E

Effective affordable mixing

Features

- Speed controlled by Scale
- Adjustable speed range from 0 to 80 rpm

Applications

Laboratory Tube Rotators are ideal for Biochemistry, Molecular Biology and Clinical applications. Provides a gentle but effective mixing of biological samples in 1.5mL to 50mL micro tubes, used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics etc.



Plastic Clips

Used with MX-RL-E/MX-RD-E Disk/Rotisserie accessory, for 50mL, 15mL and 1.5mL



Specifications	MX-RD-E
Motor type	DC motor
Speed range	0-80rpm
Speed display	Scale
Operation mode	Continuous
Voltage,Frequency	100-240V, 50/60Hz
Power	20W
Dimension [DXWXH]	300X220X310mm
Permissible ambient temperature and humidity	5-40°C,80%RH
Protection class	IP21



MX-RL-E

Economical Rotator



Features

- Provides gentle but effective rolling
- Adjustable speed range from 0 to 80 rpm

Applications

Rotisserie Horizontal Tube Rotator provides slow to vigorous mixing of your biological samples and solutions with advanced features you expect. Used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostic tests etc. Not recommended for extended use in incubators or cold rooms.

Plastic Clips

Used with MX-RL-E/MX-RD-E Disk/Rotisserie accessory, for 50mL, 15mL and 1.5mL



Specifications	MX-RL-E
Motor type	DC motor
Speed range	0-80rpm
Speed display	Scale
Operation mode	Continuous
Voltage,Frequency	100-240V, 50/60Hz
Power	20W
Dimension [DXWXH]	150X530X90mm
Permissible ambient temperature and humidity	5-40 °C,80%RH
Protection class	IP21

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: dbe@nt-rt.ru || сайт: <http://dlab.nt-rt.ru/>